



FOOT SPA MASSAGER

USER MANUAL



LC-8026

CATALOGUE

- the manual of foot health 1
- product structure 2
- usage and caution 3
- the commendatory of the function 4
- technology parameter and problem treatment 5

TECHNOLOGY PARAMETER AND PROBLEM TREATMENT

Technology parameter

Rated voltage: AC 110V	Rated power: 500W
Capability: 9L	Net weight: 3.0kg
The highest temperature: 47 °C	Frequency: 60Hz
The power line length: 2M	Specification: 51.2X31.6X41.6cm

Safeguard

The safeguard system of temperature control	Heat-resistant and aging-resistance rubber
The parts all pass the quantity safe certificate	

Problem treatment

Problem	Check	Finish
Don't heat	1. The power is off 2. There is no water when heating	1. Make sure the power is on 2. Put in water
Don't work or stop in the middletime of work	1. The temperature is too high 2. The work time is up 3. Forget to start when the power is on	1. Cold down the water 2. Restart 3. operate according to manual

THE COMMENDATORY OF THE FUNCTION

The detail of function commendatory:

- 1. A key to start**
Connect the power, the display shows OFF, press a key to start button, the machine work and all feature is turned on, you can turn off unnecessary functions selection
- 2. Dual Thermostat**
Temperature control switch, dual control temperature
- 3. Temperature regulation**
Adding water turned heated, the water temperature will be heated to around 42 degrees Celsius automatic thermostat can also be adjusted to a suitable temperature, allowing your feet to fully enjoy the hot foot bath.
- 4. Roller Massage**
A full range of health care for foot massag
- 5. Basket Design**
Look more warm and beautiful and more convenient to use
- 6. Insulation full cover**
Dustproof, improve the heating rate
- 6. Insulation full cover**
- 7. Medicine device**
Bags of herbs specific storage position, can be medicine foot bath.

The common sense of foot massager

1. Before using the machine, please put the machine on balanced level, and fix the plug, use when fixing the place.
2. In principle, it can be used two times every day, in the morning and in the evening, 30 minutes each time, the temperature should be 40 degree, and after you use it for one to two months, you will find that some symptom is eliminated.
3. After the hot wave massager, the pore of all the body is in the status of open, so do not let the electrical fan beat your body, to prevent headache.
4. When using the foot massager every time, the machine will improve the metalism, so please drink water more, in common, before and after using this machine, please drink 250ml warm water.
5. When you open the function key of gas wave and ozone, the dirt will deposit on the inner of the machine, so please wash clean after use it.
6. It will be fit to use it one hour before or one hour after the dinner. After foot massager, please do not stand and leave in time, please breath for three times, and close the eyes for 2 to 3 minutes, then walk away.

THE MANUAL OF FOOT HEALTH


Foot is the slowest part of the body, it has fifty-two bones, sixty-six arthrosis, forty muscle and two hundred ligament. People have four basis, they are ear, nose, breast, foot. And foot is one of the four important basis. Foot is the basis of the body, if it is not normal, and it will bring the abnormality of the organism. So foot is the basis of the health of the body, it is the second heart of the body.

Because of the special organization of the body, all the apparatus has nerve to foot. The medical of foot flash is the important heritage of the medical, and it has the history of 3000 years. It simulate the point that all the apparatus has in foot, then bring nerve reflect, active the function of induce apparatus, and improve the circulation of blood, adjust the maladjustment of incretion, balance the blood. It can have good curative effect for the following disease: high blood pressure, cough, asthma, waist and leg ache, rheumatism, nephropathy, tummy bug, dysmenorrhea, headache, insomnia, fever, diabetes, and so on.

The medical of foot massager can improve the blood circulation, slow the nerve of the body, excited the latent enginery, simulate the reproduction ability of cell organise, adjust the status of the imbalance of the body, make the apparatus resume itself and keep the normal function, in order to have the function of keep fit oneself and keep long life.

Specially commend to the following people:

1. the people of insomnia, nerve decline, genius nerve.
2. the people of blood circulation badness, accidie.
3. the people of leg paralytic, waist and back Ache.
4. the people of tummy, bug, digest badness.
5. the people of immunity is low, headache.
6. the people of nephric badness, sex slow.
7. the people of exercise lack, the old and healthless, arm and leg coldness.



PRODUCT STRUCTURE



- Button start
- Oxygen bubbles
- basket design
- Dual temperature control
- Temperature regulation
- Full cover keep warm
- Roller massage
- Medicine box device

USAGE AND CAUTION

USAGE

1. this foot bath not allow used in case without water
2. in order to improve the heating speed in winter, spa massager can be put hot water around 38 degrees Celsius.
3. when using the foot bath, first flat it, then power, then turn on the function
4. Press a key to start, all the features turned on, the water temperature begins to rise, it will be automatic thermostat when the water temperature reaches about 42 °C, then you can start use.
5. click bond bubble, bubbles off, and then click the function key, the bubble function is turned on.
6. press the red function key, the red light massage function turn off, press the function key, the red light is turned on.
7. After the foot bath is complete, click the start button, turn off all functions, unplug the power cord, drained the water inside, Wipe with a cloth, and the machine is stored in a dry place.

CAUTION

- Don't use the massager beyond 40 minutes for one time.
- Stop using when you don't feel good.
- To prevent the water to make the controller wet.
- It has to pull out the plug when there is a problem in draining, cleaning and using.
- Please pull out the plug immediately when you find a problem of the massager. Don't repair it by yourself. You have to take the machine to the professional repairman.
- Please put the massager in the dry condition if you don't use the machine in a long time.
- Don't put the massager in the water totally. And don't use like alcohol, corrosive chemical material, acetone to clean the machine.
- Pregnant, high blood pressure sufferer, heart disease sufferer and another serious disease sufferer use this massager seriously